

AdvanTAge™ PCR Cloning Kit Protocol-at-a-Glance (PT3067-2)

Please read the *User Manual* before using this abbreviated protocol.

PCR Amplification

- Generate PCR products using your own protocol. Be sure to finish with a final 7–10 min extension step.
- 2. Analyze 10 µl of each PCR sample by agarose/ EtBr gel electrophoresis.
- 3. Quantify the amount of DNA for each product.

Cloning into pT-Adv

1. For each desired ligation reaction, mix: Sterile $\rm H_2O$ 5 $\rm \mu I$ 10X Ligation buffer 1 $\rm \mu I$ pT-Adv (25 $\rm ng/\mu I$) 2 $\rm \mu I$ Fresh PCR product (~10 $\rm ng$) 1 $\rm \mu I$ Total 10 $\rm \mu I$

- 2. Incubate ligation reaction mixture at 14°C for at least 4 hr (preferably overnight).
- 3. Centrifuge ligation reaction mixture briefly and place on ice.

Transformation

- 1. Thaw appropriate number of tubes of TOP10F' *E. coli* competent cells on ice.
- Pipet 1–2 μl of each ligation reaction into the competent cells. Mix by gently tapping.
 (PR01083)

CLONTECH Laboratories, Inc.

Tel: 800-662-2566 (CLON) • 650-424-8222 Internet: www.clontech.com
Technical Support Fax: 650-424-1064 E-mail: tech@clontech.com

AdvanTAge™ PCR Cloning Kit Protocol-at-a-Glance, continued

Transformation continued

- 3. Incubate the tubes on ice for 30 min.
- Heat shock for exactly 30 sec in a 42°C water bath. Do not mix.
- 5. Place the tubes on ice for 2 min.
- 6. Add 250 µl of SOC medium to each tube.
- 7. Incubate tubes at 37°C for exactly 1 hr with shaking (225 rpm).
- 8. Place the tubes on ice.

Analysis

- Plate 50 μl and 200 μl from each transformation on separate LB/X-Gal/IPTG plates containing 50 μg/ml of either ampicillin or kanamycin.
- 2. Incubate at 37°C for at least 18 hr.
- 3. Shift plates to 4°C for 2–3 hr for color development.
- 4. Analyze white transformants for the presence and orientation of inserts by restriction mapping or sequencing.

The PCR process is covered by patents owned by Hoffmann-La Roche, Inc., and F. Hoffmann-La Roche, A. G.

Products covered by U.S. Patents #5,487,993 & #5,827,657 and European Patent #0550693.

© 2000, CLONTECH Laboratories, Inc.

Protocol # PT3067-2 Version # PR01083 TEL: 800-662-2566 (CLON) or 650-424-8222

page 2