

## M13 Putative Plasmid Screen PCR Worksheet (50 µl rxn volumes)

Sample(s): \_\_\_\_\_ Name: \_\_\_\_\_ Date: \_\_\_\_\_

Cocktail mixture for n = \_\_\_\_\_ (Rem: Add one more for negative control)

Component	Volume/rxn (µl)	Times n =	Total Volume Added
10X PCR Buffer	5.0		
dNTPs (10 mM)	4.0		
M13F (50 µM)	0.5		
M13R (50 µM)	0.5		
BSA (10 mg/ml)	1.0		
Mg <sup>++</sup> (25 mM)	5.0		
Taq 95 (5 U/µl)	1.0		
Subtotal =	17		
Optima Water	(50 - subtotal) = 33		
Total cocktail vol = subtotal + H <sub>2</sub> O =			
Cocktail per tube = 50 =			

Mix cocktail including Opt. H<sub>2</sub>O first, aliquot to PCR tubes, then add plasmid template from isolated colony by touching with toothpick and then rubbing tip on tube wall just submerged in PCR mix.

Method Name: \_\_\_\_\_ (Re: on ABI 9700)

Comments: