

BIOLOGY 346

MICROBIOLOGY LAB MANUAL

CONTENTS

Laboratory Safety	1-2
Microscopy	Handout
Ocular and Stage Micrometers	3-5
Simple Staining and Bacterial Morphology	6-8
Hanging Drop Technique	9-10
Differential Staining	11-16
Dilutions and Aseptic Technique	17-21
Media Preparation	22-26
Isolation of Pure Cultures	27-32
Bacterial Enumeration	33-40
Environmental Influences on Microbial Growth	41-46
Isolation of <i>Streptomyces</i> from Soil	47-49
Selective and Differential Media	50-52
Microbes of the Body: Gram-positive Cocci	53-57
Microbes of the Body: The Enterics	58-64
Multitest Identification Kits	65
Bacterial Genetics	66-68
Bactericidal Effects of Ultraviolet Radiation	69-70

Ames Test	71-73
Control of Microbial Growth: Disinfectants and Antiseptics	74-76
Control of Microbial Growth: Antibiotics	77-81
Bacterial Examination of Water	82-88
Food Microbiology	89-90
Snyder Test	91-92
Chromogenic Substrate Test	93-96

LABORATORY SAFETY

Safety in a microbiology laboratory is very important, as some of the organisms being studied are **potential** pathogens. This means that although they may not cause disease in a normal healthy host, they might if the host were compromised. A human host can be compromised in a number of different ways: cuts and wounds; lowered resistance due to disease, surgery or stress (including the stress of lecture exams); and immune-system disability.

In addition to biohazardous material, there are some chemicals used in this laboratory which are potentially harmful. Finally, many procedures will involve glassware, open flame and sharp objects which can cause damage if used improperly.

To avoid the problems that might occur, the following precautions should be taken for the safety and convenience of everyone working in the laboratory:

1. Strict attention should be given to all instructions, and any precautionary statements given in laboratory protocols should be followed.
2. The location, purpose and use of emergency safety equipment should be known before starting any assignments.
3. Do not eat, drink, smoke, or chew pencils or pens in the laboratory.
4. Laboratory bench-tops should be wiped with a disinfectant before beginning and after finishing laboratory work.
5. Hands should be washed after any laboratory work.
6. Long hair should be tied back so that it does not catch fire in the Bunsen burner flame and does not fall into sterile media.
7. Gas burners must be turned off at the end of the laboratory period.

- 8.** All cultures should be properly labeled and properly discarded after use. No slides, tubes, culture plates, flasks, or other containers with viable organisms should be washed or discarded until after they have been properly sterilized.
- 9.** Spills of chemical reagents or material containing viable organisms should be reported to the laboratory instructors and cleaned immediately.

OCULAR and STAGE MICROMETERS

Size is one of the most important physical features employed in the identification and characterization of an organism. The exact size of a microorganism can only be determined by utilizing a calibrated **ocular micrometer**. An ocular micrometer (Fig. 2-1) is a glass disc on which a series of uniformly spaced lines has been inscribed. The ocular micrometer is placed in one of the eyepieces of the microscope; however, the distance between the etched lines depends upon the objective lens used to view the specimen. In order to determine the precise distance between the lines of an ocular micrometer, it must be calibrated with a **stage micrometer** (Fig. 2-2). The inscribed lines on a stage micrometer are exactly 0.01 mm (or 10 μm) apart.

In order to calibrate the ocular micrometer for a particular objective lens, the ocular and stage micrometers are superimposed, and the number of ocular graduations per stage micrometer graduation is determined. Figure 2-3 indicates that six ocular micrometer graduations fit between two stage micrometer graduations; therefore, one space of the ocular micrometer is equal to $10\ \mu\text{m}/6$ or $1.66\ \mu\text{m}$. In order to determine the dimensions of an organism, the number of graduations occupied by the organism (Fig. 2-4) is counted and multiplied by the distance between graduations. For example, if an organism occupied the space of seven graduations, this particular dimension would be 7×1.66 or $11.6\ \mu\text{m}$.

Very seldom will an exact number of ocular micrometer lines fit between two stage micrometer lines. If that is the case, the number of stage micrometer lines is divided by the total number of ocular micrometer lines in order to determine one ocular graduation. **EXAMPLE:** 25 graduations on the ocular micrometer precisely match 4 graduations on the stage micrometer. Remembering that the graduations of the stage micrometer are $10\ \mu\text{m}$ apart, one ocular graduation = $40/25 = 1.6\ \mu\text{m}$.




Figure 2-1. Ocular micrometer.




Figure 2-3. Ocular micrometer superimposed on stage micrometer.




Figure 2-2. Stage micrometer.




Figure 2-4. Use of the ocular micrometer to measure organisms.

Material:

1. Ocular micrometer
2. Stage micrometer
3. Prepared slides of bacteria

Procedure:

1. Obtain a stage micrometer from the instructor. Place the stage micrometer on the stage of the microscope.
2. Rotate one of the microscope eyepieces until the lines of the ocular micrometer are parallel with those of the stage micrometer.

3. Match the lines on the left edges of the two micrometers by moving the stage micrometer so that the graduations of the ocular micrometer are superimposed over those of the stage micrometer.
4. Determine the number of ocular micrometer spaces that fall within a given number of stage micrometer spaces.
5. Calculate the distance between each ocular graduation by using the following formula:

$$1 \text{ ocular micrometer space } (\mu\text{m}) = \frac{\text{x spaces on the stage micrometer}}{\text{y spaces on the ocular micrometer}}$$

6. Repeat the procedure for the 10X and 40X objectives and record results. **Estimate** the calibration of the 100X oil-immersion objective lens by dividing the value for the 10X objective by 10.

Results:

One ocular micrometer
graduation equal to

4X objective	
10X objective	
40X objective	
100X objective	

SIMPLE STAINING and BACTERIAL MORPHOLOGY



Since bacteria are typically 0.2 to 2.0 μm on longest dimension, they are difficult to see with the light microscope. Bacterial cells in the "living state" were first observed by Leeuwenhoek in drops of fluid. Wet mounts were commonly used by Pasteur, Lister and many of the early workers in the field of microbiology. The method was adequate and even desirable for certain observations, but the constant motion of cells, due either to natural motility or Brownian movement, made accurate morphological studies very difficult. Robert Koch, during his classical research on anthrax, prepared thin films of bacteria on glass and allowed the smears to dry. Although the anthrax bacilli were motionless and unchanged in general shape, the more detailed structure of the cells was still obscure, due to their colorless and transparent nature. The simplest way to increase contrast with the surrounding media was to use dyes which are taken up by the bacterial cells. In 1875, Weigert, a German contemporary of Koch, found that the dye, methyl violet, would color or "stain" bacterial cells in some tissue preparations. This use of dyes to render bacterial cells easily visible under the microscope was adopted by Koch and his many students, and soon became one of the most widely used fundamental techniques of microbiology.

Prior to staining with dyes, cells must be "fixed." Fixation is a process of killing, immobilizing and preserving the bacterial cell. For bacteria, heat fixation is most common; however, chemicals such as formaldehyde, glutaraldehyde, acids and alcohols can be used.

Dyes used to stain bacterial cells are organic compounds which have affinity for specific cellular components. Many dyes are positively charged (cationic) and combine strongly with negatively charged cellular materials such as nucleic acids and acidic polysaccharides. Other dyes are anionic (negatively charged) and combine with positively charged constituents such as proteins. A dye is composed of two components:

1. An **auxochrome** group which in itself does not produce color but gives the dye its acidic or basic properties. Examples are NH_2 ; OH ; OCH_3 ; I ; Br ; and Cl .

2. A **chromophore** group which imparts color to the dye molecule.
Examples are



Simple staining is merely the use of a dye to increase the contrast of cells for microscopy. As an example, a simple stain would be used to detect the presence of bacteria in some natural material such as urine or water. It must be remembered that a simple stain alone is **not** useful as an identification tool. The purpose of this exercise is to illustrate the use of simple stains in the study of bacterial morphology.

Material:

1. Bacterial cultures
2. Crystal violet, safranin

Procedure:

1. Clean microscope slides with soap and water, rinse, and then blot dry.
2. Using an inoculating loop, put a drop of water on the slide.

3. Flame, i.e. sterilize, an inoculating loop and using **aseptic technique**, remove a small amount of bacteria from some solid medium and mix with the drop of water. The smear should be about the size of a dime and must be fairly dilute. If the smear is too dense, the morphology of individual cells will be impossible to determine.
4. Allow the smear to **air-dry**.
5. Heat-fix the slide by passing it through the flame of a Bunsen burner two or three times, or until the slide is slightly warm when touched to the back of the hand. Heat coagulates the protein of the cells so that the cells stick to the glass surface and are not washed off during the staining and rinsing procedures. Do not overheat the slide, however, as this may cause distortion of cell shape and uneven staining.
6. Allow the slide to cool, then flood with a stain - either safranin for 2 min or crystal violet for 1 min.
7. At the end of the staining time, pour off the stain and wash the slide gently but thoroughly with tap water. Blot dry with a paper towel.
8. Examine the dry smear using low-power, high-power and oil-immersion objectives. Describe the cellular morphology of each bacterial species.

HANGING DROP TECHNIQUE

Most bacterial microscopic preparations kill the organisms. The hanging drop technique allows you to observe the size, arrangement, and shape of **living** cells and to determine motility.

A thin ring of petroleum jelly is applied to the four edges on one side of a cover slip. A drop of pond scum or a drop from a bacterial broth culture is then placed in the center of the cover slip. A concave (depression) microscope slide is carefully placed over the cover slip in such a way that the drop is suspended and is undisturbed. The petroleum jelly causes the cover slip to stick to the slide. The slide preparation may now be inverted and placed under the microscope for examination. Since no stain is used and most cells are transparent, viewing is best done with as little illumination as possible. The petroleum jelly form an air-tight seal that prevents drying of the drop allowing a long period to observe cell size and shape, binary fission, and motility.

If the technique is used to determine motility, the observer must be careful to distinguish between true motility and Brownian motion. Brownian motion is caused by water molecules colliding with the organism and moving it around in an irregular jerky pattern. Organisms will appear to vibrate in place. With true motility, cells will exhibit independent movement in some consistent direction over greater distances.


FIRST PERIOD

Material:


1. Depression slide and cover slip
2. Petroleum jelly
3. Toothpick
4. Pond scum or 24-h bacterial broth culture

Procedure:


1. *Apply the petroleum jelly.* Use a toothpick to place a thin ring of petroleum jelly around the well of a depression slide. An easy short cut is to only spot the petroleum jelly on the four corners of the cover slip, though drying can then occur.



2. *Transfer the Organisms.* Use your loop to place a drop of pond scum on a cover slip. Do not spread out the drop.



3. *Invert the Slide over the cover slip.* Carefully invert the depression slide over the cover slip so the drop is centered in its well. Gently press until the petroleum jelly has created a seal between the slide and the cover slip.



4. *Observe the slide.* Correctly done, the slide should look like this from the side. Place the slide with the cover slip up on the microscope stage and observe under high dry or oil immersion.

